

# Anti-Glutaminase C Rabbit antibody



Catalog Number: 168897

## Product name

Rabbit IgG

Anti-Glutaminase C Rabbit antibody

## Applications

## Specificity

WB, ICC/IF, FC, IP

Human, Mouse, Rat

## Dilutions

## Antibody description

1:1000

Glutaminase C Rabbit monoclonal antibody

1:50

## Preparation

1:50

Antigen: Recombinant protein of human  
Glutaminase C

1:20

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M Nacl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of Glutaminase in K562,3T3,Hela cell lysates using Glutaminase Rabbit mAb(1:1000 diluted).Predicted band size:74kDa.Observed band size:68kDa.