

CBFA2T3 antibody



Catalog Number: 108897

Product name

CBFA2T3 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

CBFA2T3 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 62-65kd

Preparation

This antibody was obtained by immunization of CBFA2T3 recombinant protein (Accession Number: BC062624). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:108897(CBFA2T3 antibody) at dilution of 1:1000