

NKG2A antibody

Catalog Number: 113197

Product name

NKG2A antibody

Specificity

Human; other species not tested.

Antibody description

NKG2A Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HEK-293 cells, human liver tissue. Positive IHC detected in human lung cancer tissue. Observed molecular weight by Western-blot: 33kd

Preparation

This antibody was obtained by immunization of NKG2A recombinant protein (Accession Number: NM_007328). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113197(KLRC1 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:113197(KLRC1 antibody) at dilution of 1:100 (under 10x lens)