

# SERPINA12 antibody



Catalog Number: 115197

## Product name

SERPINA12 antibody

## Specificity

Human; other species not tested.

## Antibody description

SERPINA12 Rabbit Polyclonal antibody. Positive WB detected in human adipose tissue. Observed molecular weight by Western-blot: 47 kDa

## Preparation

This antibody was obtained by immunization of SERPINA12 recombinant protein (Accession Number: XM\_017020990). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human adipose tissue were subjected to SDS PAGE followed by western blot with Catalog No:115197(SERPINA12 Antibody) at dilution of 1:1000