

Anti-TCF-4/TCF7L2 Rabbit antibody



Catalog Number: 168597

Product name

Anti-TCF-4/TCF7L2 Rabbit antibody

Applications

WB, IHC-P, ICC/IF, FC, IP

Specificity

Human, Mouse, Rat

Dilutions

1:5000

Antibody description

TCF-4/TCF7L2 Rabbit monoclonal antibody

1:50

1:50

Preparation

Antigen: A synthetic peptide of human TCF-4

1:20

1:20

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Validations

Storage


Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG



Western blot detection of TCF-4 in K562, HeLa cell lysates using TCF-4 Rabbit mAb(1:1000 diluted). Predicted band size: 68kDa. Observed band size: 68kDa.