

# MYOZ1 antibody

Catalog Number: 112997

## Product name

MYOZ1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MYOZ1 Rabbit Polyclonal antibody. Positive IP detected in mouse skeletal muscle tissue. Positive WB detected in mouse skeletal muscle tissue, human heart tissue. Observed molecular weight by Western-blot: 32 kDa

## Preparation

This antibody was obtained by immunization of MYOZ1 recombinant protein (Accession Number: NM\_021245). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:112997(MYOZ1 Antibody) at dilution of 1:600



IP Result of anti-MYOZ1 (IP:Catalog No:112997, 4ug; Detection:Catalog No:112997 1:500) with mouse skeletal muscle tissue lysate 4000ug.