

Anti-Phospho-GAB1 (Tyr627) Rabbit antibody



Catalog Number: 166997

Product name

Anti-Phospho-GAB1 (Tyr627) Rabbit antibody

Specificity

Human, Mouse

Antibody description

Phospho-GAB1 (Tyr627) Rabbit polyclonal antibody

Preparation

Antigen: Peptide sequence around phosphorylation site of tyrosine 627 (V-E-Y(p)-L-D) derived from Human Gab1.

Formulation

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB, ICC/IF


Dilutions

WB:1/500-1/1000 ICC:1/100-1/200


Validations



Western blot detection of Gab1(Phospho-Tyr627) in A549 cells untreated or treated with EGF using Gab1(Phospho-Tyr627) Rabbit pAb (dilution 1:500, upper) or Gab1 Rabbit pAb (310130, dilution 1:1000, lower). Predicted band size: 110kDa. Observed band size: 110kDa.



Immunofluorescence staining of methanol-fixed HeLa cells using Gab1(Phospho-Tyr627) Antibody .



Western blot analysis of extracts from 293 cells untreated or treated with EGF using Gab1 (Phospho-Tyr627) Antibody .