

ZNF496 antibody



Catalog Number: 116997

Product name

ZNF496 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZNF496 Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 67 kDa

Preparation

This antibody was obtained by immunization of ZNF496 recombinant protein (Accession Number: NM_032752). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116997(ZNF496 antibody) at dilution of 1:1000