

Anti-AKAP 95 Rabbit antibody



Catalog Number: 167997

Product name

Anti-AKAP 95 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

AKAP 95 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human AKAP 95

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/1000

IHC: 1/50

ICC/IF: 1/50

Validations



Western blot detection of AKAP 95 in HeLa, A549, HL-60, U251 cell lysates using AKAP 95 Rabbit mAb(1:1000 diluted). Predicted band size: 76kDa. Observed band size: 95kDa.