

Anti-P2Y6 Rabbit antibody



Catalog Number: 168997

Product name

Anti-P2Y6 Rabbit antibody

Specificity

Human

Antibody description

P2Y6 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human P2Y6

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP

Dilutions


1:2000-1:10000

1:20-1:100

1:20-1:50

1:20

Validations



Western blot detection of P2Y6 in K562, HeLa cell lysates using P2Y6 Rabbit mAb(1:1000 diluted). Predicted band size: 36kDa. Observed band size: 36kDa.