

# HRP-conjugated GFP tag antibody



Catalog Number: 107297

## Product name

Monoclonal

HRP-conjugated GFP tag antibody

## Ig Type

## Specificity

Mouse IgG2a

Tag

## Applications

## Antibody description

ELISA, WB

HRP-conjugated GFP tag Mouse Monoclonal antibody. Positive WB detected in Recombinant protein. Observed molecular weight by Western-blot: -

## Dilutions


Recommended Dilution:

WB: 1:5000-1:50000

## Preparation

This antibody was obtained by immunization of Recombinant Protein (Accession Number: AVU73901). Purification method: Protein A Purified.

## Validations



## Formulation

PBS with 0.02% Thimerosal, 50% Glycerol, pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Recombinant protein were subjected to SDS PAGE followed by western blot with HRP-66002(GFP tag Antibody) at dilution of 1:10000