

HRP-conjugated GFP tag antibody



Catalog Number: 107297

Product name

HRP-conjugated GFP tag antibody

Specificity

Tag

Antibody description

HRP-conjugated GFP tag Mouse Monoclonal antibody. Positive WB detected in Recombinant protein. Observed molecular weight by Western-blot: -

Preparation

This antibody was obtained by immunization of Recombinant Protein (Accession Number: AVU73901). Purification method: Protein A Purified.

Formulation

PBS with 0.02% Thimerosal, 50% Glycerol, pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:5000-1:50000

Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with HRP-66002(GFP tag Antibody) at dilution of 1:10000