

SLC16A3 Polyclonal Antibody



Catalog Number: 162697

Product name

SLC16A3 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to SLC16A3

Preparation

Antigen: A synthetic peptide of human MCT4

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:200 - 1:500

IHC 1:50 - 1:100

Validations



Western blot - SLC16A3 Polyclonal Antibody

Western blot analysis of extracts of NCI-H460 cells, using SLC16A3 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.