

# Mouse S100A5 (His Tag) recombinant protein



Catalog Number: 502697

## General Information

### Gene Name Synonym

Protein S-100D; S100 calcium-binding protein A5

### Protein Construction

A DNA sequence encoding the mouse S100A5 (P63084) (Met 1-Lys 93) was expressed, with a polyhistidine tag at the N-terminus.

### Organism

Mouse

### Expression Host

E. coli

## QC Testing

### Purity

> 97 % as determined by SDS-PAGE

### Endotoxin

Please contact us for more information.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Met

### Molecular Mass

The recombinant mouse S100A5 consisting of 104 amino acids and has a calculated molecular mass

of 12.3 kDa. rmS100A5 migrates as an approximately 13 kDa band in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile 50mM Tris, 20% glycerol, pH 7.5

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Mouse S100A5 Protein (His Tag) SDS-PAGE