

# Anti-FGFR2 antibody

Catalog Number: 176697

## Product name

Anti-FGFR2 antibody

## Specificity

Human

## Antibody description

Rabbit polyclonal antibody to FGFR2

## Preparation

This antigen of this antibody was recombinant protein within human fgfr2/cd332 aa 600-750 / 821.

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, FC

## Dilutions

WB: 1:1,000

IHC-P: 1:200

## Validations




Fig1: Western blot analysis on MCF-7 cell lysates using anti- CD332 rabbit polyclonal antibodies.




Fig2: Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using anti- CD332 rabbit polyclonal antibody.