

UTP3 antibody



Catalog Number: 116697

Product name

UTP3 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

UTP3 Rabbit Polyclonal antibody. Positive WB detected in A549 cells, HeLa cells, human brain tissue, mouse thymus tissue. Positive IF detected in HepG2. Observed molecular weight by Western-blot: 70kd

Preparation

This antibody was obtained by immunization of UTP3 recombinant protein (Accession Number: NM_020368). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:20-1:200

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:116697(UTP3 antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells, using UTP3 antibody Catalog No:116697 at 1:50dilution and Rhodamine-labeled goat anti-rabbit IgG (red).