

Anti-Hck Rabbit antibody



Catalog Number: 167697

Product name

Anti-Hck Rabbit antibody

Specificity

Human

Antibody description

Hck Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Hck

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/2000-1/10000

ICC/IF: 1/50

FC: 1/20

IP: 1/20

Validations



Western blot detection of Hck in HeLa, A549, HL-60, U2OS cell lysates using Hck Rabbit mAb(1:1000 diluted). Predicted band size: 60kDa. Observed band size: 60kDa.