

# MIER2 Polyclonal Antibody



Catalog Number: 164397

## Product name

MIER2 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to MIER2

## Preparation

Antigen: A synthetic Peptide of human MIER2

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, FC

## Dilutions

WB 1:500 - 1:1000

FC 1:20 - 1:50

## Validations



Western blot - MIER2 Polyclonal Antibody

Western blot analysis of extracts of NCI-H460 cells, using MIER2 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.