

IL-29 antibody

Catalog Number: 107397

Product name

IL-29 antibody

Specificity

Human; other species not tested.

Antibody description

IL-29 Mouse Monoclonal antibody. Positive WB detected in Recombinant protein, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 22 kDa

Preparation

This antibody was obtained by immunization of IL-29 recombinant protein (Accession Number: NM_172140). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:5000-1:50000

Validations



GST tagged recombinant IL29 protein were subjected to SDS PAGE followed by western blot with Catalog No:107397(IL29 Antibody) at dilution of 1:20000



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107397(IL29 Antibody) at dilution of 1:5000