

Anti-Ferritin Light Chain (1B4) Mouse antibody



Catalog Number: 168397

Product name

Monoclonal

Anti-Ferritin Light Chain (1B4) Mouse antibody

Ig Type

Mouse IgG2a

Specificity

Human

Applications

Antibody description


ELISA

Ferritin Light Chain (1B4) Mouse monoclonal antibody

Validations

Preparation

Antigen: Purified recombinant human TAK1 protein fragments expressed in E.coli.



Formulation

PBS(pH 7.4) containing with 0.03% Proclin300.

Storage

Store at 4°C short term.Avoid freeze / thaw cycle.

Clonality

Observed Ferritin Light Chain levels in Human serum samples at different dilution.