

# Syntaxin 3 antibody

Catalog Number: 115797

## Product name

Syntaxin 3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Syntaxin 3 Rabbit Polyclonal antibody. Positive IP detected in SH-SY5Y cells. Positive WB detected in SH-SY5Y cells, human blood tissue, Neuro-2a cells. Observed molecular weight by Western-blot: 33 kDa

## Preparation

This antibody was obtained by immunization of Syntaxin 3 recombinant protein (Accession Number: XM\_005274195). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:115797 (STX3 Antibody) at dilution of 1:1000



IP Result of anti-STX3 (IP:Catalog No:115797, 4ug; Detection:Catalog No:115797 1:1000) with SH-SY5Y cells lysate 3600ug.