

Anti-ATG10 Rabbit antibody



Catalog Number: 167797

Product name

Anti-ATG10 Rabbit antibody

Specificity

Human

Antibody description

ATG10 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human ATG10

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 1/1000

IP: 1/20

Validations



Western blot detection of ATG10 in HeLa,A549,U251 cell lysates using ATG10 Rabbit mAb(1:1000 diluted).Predicted band size:25kDa.Observed band size:25kDa.