

# C9orf11 antibody



Catalog Number: 108797

## Product name

C9orf11 antibody

## Immunogen

[Human C9orf11 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

C9orf11 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in HepG2 cells, HeLa cells. Observed molecular weight by Western-blot: 27-29 kDa

## Preparation

This antibody was obtained by immunization of C9orf11 recombinant protein (Accession Number: NM\_020641). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:108797(C9orf11 antibody) at dilution of 1:2000



Immunofluorescent analysis of HepG2 cells using Catalog No:108797(C9orf11 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG