

Mouse PLAUR/CD87 (His & Fc Tag) recombinant protein



Catalog Number: 500129

General Information

Gene Name Synonym:

Protein Construction:

A DNA sequence encoding the extracellular domain (Met 1-Thr 297) of mouse PLAUR (NP_035243.1) precursor was fused with the C-terminal polyhistidine-tagged Fc region of human IgG1 at the C-terminus.

Source: Mouse

Expression Host: Human Cells

QC Testing

Activity:

Measured by its binding ability in a functional ELISA. Immobilized human uPA at 5 µg/ml (100 µl/well) can bind mouse PLAUR with a linear ranger of 1.6-40 ng/ml.

Purity:

> 97 % as determined by SDS-PAGE

Endotoxin:

< 1.0 EU per µg of the protein as determined by the LAL method

Stability:

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal:

Leu 24

Molecular Mass:

The recombinant mouse PLAUR/Fc is a disulfide-

linked homodimer after removal of the signal peptide. The reduced monomer consists of 522 amino acids and has a predicted molecular mass of 58 kDa. In SDS-PAGE under reducing conditions, the apparent molecular mass of rmPLAUR/Fc monomer is approximately 80-90 kDa due to glycosylation.

Formulation:

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage:

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution:

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE:



Mouse PLAUR / CD87 / uPAR Protein (His & Fc Tag) SDS-PAGE