

Rat RAB7A / Rab-7a (His Tag) recombinant protein



Catalog Number: 500140

General Information

Gene Name Synonym

Ras-related protein BRL-RAS; Ras-related protein p23

Protein Construction

A DNA sequence encoding the rat Rab7a (NP_076440.1) (Met1-Cys207) was expressed with a polyhistidine tag at the N-terminus.

Organism

Rat

Expression Host

E. coli

QC Testing

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant rat Rab7a consists 225 amino acids and predicts a molecular mass of 25.7 kDa.

Formulation

Lyophilized from sterile PBS, 10 % glycerol.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C.

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat RAB7A / Rab-7a Protein (His Tag)