

Rat UNC5B / UNC5H2 (Fc Tag) recombinant protein



Catalog Number: 500142

General Information

Gene Name Synonym

Protein unc-5 homolog 2; Protein unc-5 homolog B

Protein Construction

A DNA sequence encoding the rat UNC5B (NP_071543.1) (Met1-Asp373) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gly 27

Molecular Mass

The recombinant rat UNC5B/Fc is a disulfide-

linked homodimer. The reduced monomer comprises 588 amino acids and has a predicted molecular mass of 65.7 kDa. The apparent molecular mass of the protein is approximately 72-90 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat UNC5B / UNC5H2 Protein (Fc Tag) SDS-PAGE