

Rat B7-H4/B7S1/B7x (His Tag) recombinant protein



Catalog Number: 500350

General Information

Protein Construction

A DNA sequence encoding the rat Vtn1 (NP_001019415.1) (Phe29-Ser256) was expressed with a polyhistidine tag at the C-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Phe 29

Molecular Mass

The recombinant rat Vtn1 consists 239 amino

acids and predicts a molecular mass of 26.6 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C.

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Vtn1 protein SDS-PAGE