

Human Contactin 4/CNTN4 (His Tag) recombinant protein



Catalog Number: 500469

General Information

Gene Name Synonym

Brain-derived immunoglobulin superfamily protein 2

Protein Construction

A DNA sequence encoding the mature form of human CNTN4 isoform 1 (Q8I WV2-1) (Met 1-Ser 1000) was fused with a polyhistidine tag at the C-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 19

Molecular Mass

The secreted recombinant human CNTN4 consists of 993 amino acids and has a calculated molecular mass of 110 kDa. It migrates as an approximately 120-130 kDa band as predicted in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20mM Tris, 500mM NaCl, pH 8.5, 10% gly

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human Contactin 4 / CNTN4 Protein (His Tag)
SDS-PAGE