

# Human DCX (aa 45-150, GST Tag) recombinant protein



Catalog Number: 500541

## General Information

### Gene Name Synonym

Doublin; Lissencephalin-X

### Protein Construction

A DNA sequence encoding the human DCX (O43602-2) N-terminal fragment (Ala 45-Val 150) was fused with the GST tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

> 82 % as determined by SDS-PAGE

### Endotoxin

Please contact us for more information.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Met

### Molecular Mass

The recombinant human DCX (aa 45-150)/GST chimera consists of 340 amino acids and has a

predicted molecular mass of 39.4 kDa. It migrates as an approximately 36 kDa band in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile 20mM Tris, 1mM DTT, 10% glycerol, pH 7.5

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human DCX Protein (aa 45-150, GST Tag) SDS-PAGE