

# Human ALK-7 / ACVR1C (Fc Tag) recombinant protein



Catalog Number: 500577

## General Information

### Gene Name Synonym

Activin receptor type IC; Activin receptor-like kinase 7

### Protein Construction

A DNA sequence encoding the human ACVR1C (NP\_660302.2) (Met1-Glu113) was expressed with the Fc region of human IgG1 at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Leu 22

### Molecular Mass

The recombinant human ACVR1C consists 330 amino acids and predicts a molecular mass of 36.6 kDa.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



ACVR1C protein SDS-PAGE