

# Human ABHD10 (aa 53-306, His Tag) recombinant protein



Catalog Number: 500606

## General Information

### Gene Name Synonym

Alpha/beta hydrolase domain-containing protein 10

### Protein Construction

A DNA sequence encoding the human ABHD10 (Q9NUJ1)(Thr53-Asn306) was fused with a polyhistidine tag at the N-terminus.

### Organism

Human

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Purity

> 90 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

His

### Molecular Mass

The recombinant human ABHD10 consists of 271

amino acids and has a calculated molecular mass of 30.3 kDa. The recombinant protein migrates as an approximately 32 kDa band in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile 20mM Tris, 500mM NaCl, 10% glycerol, pH 8.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide

### Storage


Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human ABHD10 Protein (aa 53-306, His Tag)  
SDS-PAGE