

Human CD44 (His Tag) recombinant protein



Catalog Number: 500618

General Information

Gene Name Synonym

CDw44; Epican; Extracellular matrix receptor III; GP90 lymphocyte homing/adhesion receptor; HUTCH-I; Heparan sulfate proteoglycan; Hermes antigen; Hyaluronate receptor; Phagocytic glycoprotein 1; Phagocytic glycoprotein I

Protein Construction

A DNA sequence encoding the N-terminal fragment (Met 1-Pro 220) of human CD44 (P16070-1) extracellular domain was fused with a polyhistidine tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 21

Molecular Mass

The recombinant human CD44 consists of 211 amino acids and has a predicted molecular mass of 23.4 kDa. As a result of glycosylation, the apparent molecular mass of rhCD44 is approximately 43-48 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human CD44 Protein (His Tag) SDS-PAGE