

H12N3 HA (His Tag) recombinant protein



Catalog Number: 500644

General Information

Protein Construction

A DNA sequence encoding the Influenza A virus (A/bar headed goose/Mongolia/143/2005(H12N3)) hemagglutinin (ACV86810.1), (Met1-Arg342) termed as HA1, was expressed with a polyhistidine tag at the C-terminus.

Organism

H12N3

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Tyr 17

Molecular Mass

The recombinant HA1 subunit of the Influenza A virus hemagglutinin (A/bar headed goose/Mongolia/143/2005(H12N3)) consists 337

amino acids and predicts a molecular mass of 38.4 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H12N3 (A/bar headed goose/Mongolia/143/2005) Hemagglutinin / HA1 Protein (His Tag)