

Rat OLR1/LOX1 (Fc Tag) recombinant protein



Catalog Number: 500692

General Information

Gene Name Synonym

Lectin-like oxidized LDL receptor 1; Lectin-type oxidized LDL receptor 1; Oxidized low-density lipoprotein receptor 1, soluble form

Protein Construction

A DNA sequence encoding the rat OLR1 (O70156) (Leu60-Gln364) was expressed with Fc region of human IgG1 at the N-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Purity

(90.4+8.5) % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu

Molecular Mass

The recombinant rat OLR1 comprises 565 amino acids and predicts a molecular mass of 63.9 kDa. The apparent molecular mass of the recombinant protein is approximately 66 and 35 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat OLR1 / LOX1 Protein (Fc Tag) SDS-PAGE