

# Human HBP1 (GST Tag) recombinant protein



Catalog Number: 500761

## General Information

### Gene Name Synonym

HMG box transcription factor 1; High mobility group box transcription factor 1

### Protein Construction

A DNA sequence encoding the mature form of human HBP1 (O60381-1) (Pro208-Phe345) was fused with the GST tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

Please contact us for more information.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Met

### Molecular Mass

The recombinant human HBP1 /GST chimera consists of 372 amino acids and has a predicted molecular mass of 42.5 kDa. It migrates as an approximately 43 KDa band in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human HBP1 Protein (GST Tag) SDS-PAGE