

Human COL9A1 (Fc Tag) recombinant protein



Catalog Number: 500763

General Information

Protein Construction

A DNA sequence encoding the human COL9A1 (P20849-3) (Met1-Pro328) was fused with Fc region of mouse IgG at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ala 24

Molecular Mass

The recombinant human COL9A1/mFc is a disulfide-linked homodimer. The reduced monomer comprises 539 amino acids and has a predicted molecular mass of 59.2 kDa. The

apparent molecular mass of the protein is approximately 62-72 in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human COL9A1 Protein (Fc Tag) SDS-PAGE