

Human EG-VEGF / prokineticin-1 (His Tag) recombinant protein



Catalog Number: 500929

General Information

Gene Name Synonym

Endocrine-gland-derived vascular endothelial growth factor; Mambakine

Protein Construction

A DNA sequence encoding the human EG-VEGF (Met 1-Phe105) (P58294) was expressed, with a C-terminal polyhistidine tag.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 89 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ala 20

Molecular Mass

The secreted recombinant human EG-VEGF

consists of 96 amino acids and predicts a molecular mass of 11 KDa. The apparent molecular mass of the protein is approximately 15 Kda in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile 20mM Tris, 500mM NaCl, pH 7.4, 0.02% Tween-80, 10% gly, 1mM DTT
1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human EG-VEGF / prokineticin-1 Protein (His Tag) SDS-PAGE