

H6N8 HA (His Tag) recombinant protein



Catalog Number: 501003

General Information

Protein Construction

A DNA sequence encoding the Influenza A virus (A/mallard/Ohio/217/1998(H6N8)) hemagglutinin (ABO52049.1) (Met1-Gln529) was expressed with a C-terminal polyhistidine tag.

Organism

H6N8

Expression Host

Baculovirus-Insect Cells

QC Testing

Activity

1. Measured by its ability to bind with Neu5Aca2-3Galb1-4GlcNAcb-PAA-biotin (01-077) using the Octet RED System.
2. Measured by its ability to bind with Neu5Aca2-6GalNAca-PAA-biotin (01-059) using the Octet RED System.
3. Measured by its ability to agglutinate guinea pig red blood cells. HA titer is 0.05-0.3 $\mu\text{g/mL}$ for 1% GRBC.

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 17

Molecular Mass

The recombinant hemagglutinin of Influenza A virus (A/mallard/Ohio/217/1998(H6N8)) comprises 524 amino acids and has a predicted molecular mass of 59.3 kDa. The apparent molecular mass of the protein is approximately 56.8 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10 % glycerol, pH 8.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H6N8 (A/mallard/Ohio/217/1998)
Hemagglutinin / HA Protein (His Tag) SDS-PAGE