

Mouse FGF21 (His Tag) recombinant protein



Catalog Number: 501054

General Information

Protein Construction

A DNA sequence encoding the mouse FGF21 (Q9JJN1) (Met 1-Ser 210) was expressed, with a polyhistidine tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ala 29

Molecular Mass

The secreted recombinant mouse FGF21 consists of 193 amino acids and has a calculated molecular mass of 21.3 kDa. As a result of glycosylation, the apparent molecular mass of rm FGF21 is

approximately 25 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse FGF21 / Fibroblast Growth Factor 21 Protein (His Tag) SDS-PAGE