

Human CD164 (His Tag) recombinant protein



Catalog Number: 501127

General Information

Gene Name Synonym

Endolyn; Multi-glycosylated core protein 24

Protein Construction

A DNA sequence encoding the human CD164 (Q04900-1) extracellular domain (Met 1-Asp 162) was fused with a polyhistidine tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 24

Molecular Mass

The recombinant human CD164 consists of 150

amino acids and has a predicted molecular mass of 15.9 kDa. As a result of different glycosylation, the apparent molecular mass of rhCD164 is approximately 40-110 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human CD164 / Endolyn Protein (His Tag) SDS-PAGE