

Human PPIL2 (aa 280-457, His Tag) recombinant protein



Catalog Number: 501130

General Information

Gene Name Synonym

CYC4; Cyclophilin-60; Rotamase PPIL2

Protein Construction

A DNA sequence encoding the mature form of human PPIL2 (Q13356-2) (Gly280-Pro457) was expressed with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human PPIL2 consists of 193

amino acids and predicts a molecular mass of 21.9 KDa. It migrates as an approximately 24 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human PPIL2 Protein (aa 280-457, His Tag) SDS-PAGE