

Rat BAFFR / TNFRSF13C / CD268 (His Tag) recombinant protein



Catalog Number: 501147

General Information

Protein Construction

A DNA sequence encoding the rat TNFRSF13C (D3ZNT3) (Ser10-Ala71) was expressed with a polyhistidine tag at the C-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ser 10

Molecular Mass

The recombinant rat TNFRSF13C comprises 73 amino acids and predicts a molecular mass of 8.3 kDa. The apparent molecular mass of the

recombinant protein is approximately 20 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat BAFFR / TNFRSF13C Protein (His Tag) SDS-PAGE