

Mouse LAMP2/CD107b (His Tag) recombinant protein



Catalog Number: 501194

General Information

Gene Name Synonym

CD107 antigen-like family member B; Lysosomal membrane glycoprotein type B

Protein Construction

A DNA sequence encoding the mouse LAMP2 (P17047-1) extracellular domain (Leu 26-Asn 379) was fused with a polyhistidine tag at the C-terminus and a signal peptide at the N-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 97 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Leu 26

Molecular Mass

The secreted recombinant mouse LAMP2

comprises 365 amino acids and has a calculated molecular mass of 40.6 kDa. As a result of glycosylation, the recombinant protein migrates as an approximately 70-80 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse LAMP2 / CD107b Protein (His Tag) SDS-PAGE