

Human CDC42 (GST Tag) recombinant protein



Catalog Number: 501214

General Information

Gene Name Synonym

G25K GTP-binding protein

Protein Construction

A DNA sequence encoding the mature form of human CDC42 isoform 2 (P60953-2) (Met 1-Cys 188) was fused with the GST tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human CDC42/GST chimera

consists of 423 amino acids and has a predicted molecular mass of 48.1 kDa. It migrates as an approximately 44 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20mM Tris, 0.15M NaCl, 0.5mM GSH, pH 8.0

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human CDC42 / G25K Protein (GST Tag) SDS-PAGE