

Mouse CDC2 Kinase / CDK1 (His & GST Tag) recombinant protein



Catalog Number: 501223

General Information

Gene Name Synonym:

Protein Construction:

A DNA sequence encoding the mouse CDK1 (P11440) (Met1-Met297) was expressed with the N-terminal polyhistidine-tagged GST tag at the N-terminus.

Source: Mouse

Expression Host: Baculovirus-Insect Cells

QC Testing

Activity:

Kinase activity untested

Purity:

> 85 % as determined by SDS-PAGE

Endotoxin:

< 1.0 EU per µg of the protein as determined by the LAL method

Stability:

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal:

Met

Molecular Mass:

The recombinant mouse CDK1/GST chimera consists of 534 amino acids and has a calculated

molecular mass of 61.9 kDa. The recombinant protein migrates as an approximately 57 kDa band in SDS-PAGE under reducing conditions.

Formulation:

Supplied as sterile 20mM Tris, 500mM NaCl, 10% glycerol, pH 8.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage:

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution:

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE:



Mouse CDC2 Kinase / CDK1 Protein (His & GST Tag) SDS-PAGE