

Human CDK4 (GST Tag) recombinant protein



Catalog Number: 501272

General Information

Gene Name Synonym:

Protein Construction:

A DNA sequence encoding the human CDK4 (NP_000066.1) (Met 1-Glu 303) was fused with the GST tag at the N-terminus.

Source: Human

Expression Host: Baculovirus-Insect Cells

QC Testing

Activity:

No Kinase Activity

Purity:

> 90 % as determined by SDS-PAGE

Endotoxin:

< 1.0 EU per μg of the protein as determined by the LAL method

Stability:

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal:

Met

Molecular Mass:

The recombinant human CDK4/GST chimera consists of 528 amino acids and predicts a molecular mass of 60 kDa. It migrates as an

approximately 55 kDa band in SDS-PAGE under reducing conditions.

Formulation:

Lyophilized from sterile 50mM Tris, 100mM NaCl, 10% gly, 0.5mM PMSF, pH 8.0

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage:

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution:

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE:



Human CDK4 / CMM3 Protein (GST Tag) SDS-PAGE