

H5N9 HA (His Tag) recombinant protein



Catalog Number: 501293

General Information

Protein Construction

A DNA sequence encoding the Influenza A virus (A/chicken/Italy/22A/ 1998(H5N9)) hemagglutinin (ABR37720.1) (Met1-Arg342) was expressed with a C-terminal polyhistidine tag.

Organism

H5N9

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 17

Molecular Mass

The recombinant hemagglutinin of Influenza A virus (A/chicken/Italy/22A/ 1998(H5N9)) comprises 337 amino acids and has a predicted

molecular mass of 38.1 kDa. The apparent molecular mass of the protein is approximately 48.5 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H5N9 (A/chicken/Italy/22A/1998) Hemagglutinin / HA Protein (His Tag) SDS-PAGE