

Mouse DOPA Decarboxylase/DDC (His Tag) recombinant protein



Catalog Number: 501294

General Information

Gene Name Synonym

DOPA decarboxylase

Protein Construction

A DNA sequence encoding the mouse DDC (O88533) (Met 1-Glu 480) was expressed with a C-terminal polyhistidine tag.

Organism

Mouse

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant mouse DDC consists of 490

amino acids and has a calculated molecular mass of 55.2 kDa. It migrates as an approximately 50 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 50mM Tris, 100mM NaCl, pH 8.0, 20% gly, 3mM DTT

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse DOPA Decarboxylase / DDC Protein (His Tag) SDS-PAGE