

H6N1 HA (His Tag) recombinant protein



Catalog Number: 501345

General Information

Protein Construction

A DNA sequence encoding the N-terminal segment (Met 1-Arg 344) of the influenza A hemagglutinin (A/northern shoveler/California/HKWF115/2007 (H6N1)) (ACE81692.1), termed as HA1 was expressed, fused with a C-terminal polyhistidine tag.

Organism

H6N1

Expression Host

Human Cells

QC Testing

Purity

> 92 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 17

Molecular Mass

The secreted recombinant influenza A H6N1 HA1 subunit (A/northern shoveler/California/HKWF115/2007 (H6N1)) comprises 339 amino acids and has a predicted

molecular mass of 38 kDa. As a result of glycosylation, the apparent molecular mass of the protein is approximately 45-55 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H6N1 (A/northern shoveler/California/HKWF115/2007) Hemagglutinin Protein (HA1 Subunit) (His Tag) SDS-PAGE