

# Mouse TLR2/CD282 (His Tag) recombinant protein



Catalog Number: 501365

## General Information

### Protein Construction

A DNA sequence encoding the mouse TLR2 (Q9QUN7) (Met1-Gln587) was fused with a polyhistidine tag at the C-terminus.

### Organism

Mouse

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Purity

> 90 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu$ g of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Gln 25

### Molecular Mass

The recombinant mouse TLR2 consists of 573 amino acids and has a calculated molecular mass of 64.8 kDa. The recombinant protein migrates as

an approximately 65 kDa band in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile 20mM Tris, 500mM NaCl, pH 7.4, 10% gly

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Mouse TLR2 Protein (His Tag) SDS-PAGE